

PRODUCT DATA SHEET

Description	Monoclonal antibody to Triiodothyronine
REF	10-13
Specificity	Reacts with T3
Clone	3B10, IgG1
Host Animal	Mouse
Immunogen	T3-KLH
Cross-reactivity	Thyroxin: 0.4%
Source	Ascitic fluid
Format	Purified, liquid
Purification	Chromatography on protein G sepharose; purity is tested by electrophoresis
Concentration	Lot-specific
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium azide
Long term storage	-18...-22°C
Handling	Avoid freeze/thaw cycles
Application	ELISA