

**PRODUCT DATA SHEET**

<b>Description</b>	<b>Monoclonal antibody to Testosterone</b>
<b>REF</b>	<b>10-03</b>
Specificity	Reacts with Testosterone
Clone	3F6, IgG1
Host Animal	Mouse
Source	Ascitic fluid
Format	Purified, liquid
Purification	Chromatography on Protein G Sepharose, purity is tested by electrophoresis
Concentration	Lot-specific
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Long term storage	-18...-22°C
Transportation conditions	Short term at +4°C
Handling	Avoid freeze/thaw cycles
Application	ELISA
Reference	<b>Experimental approach to study of the gonadotropin-releasing hormone function by knockdown of its expression with RNA interference in ontogenesis</b> Dygalo NN, Biological Sciences, V. 426, № 1, 2009, 213-215