

PRODUCT DATA SHEET

Description	Monoclonal antibody to CA-125
REF	10-12
Specificity	Reacts with human CA-125 (M-11-like)
Clone	2B6, IgG1
Host Animal	Mouse
Immunogen	CA-125
Cross-reactivity	Does not cross-react with CA 19-9, CA 15-3
Source	Ascitic fluid
Format	Purified, liquid
Purification	Chromatography on Protein G Sepharose, purity is tested by electrophoresis.
Concentration	Lot-specific
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium azide
Long term storage	-18...-22°C
Handling	Avoid freeze/thaw cycles
Application	ELISA, use as a solid-phase antibody together with Mab 3C8