

**PRODUCT DATA SHEET**

<b>Description</b>	<b>Ovarian tumor antigen CA-125 (Cancer antigen 125), Human</b>
<b>REF</b>	<b>11-02</b>
Source	OVCAR3 ovarian adenocarcinoma supernatants
Format	Purified, liquid
Purification	Gel-filtration chromatography after the ammonium sulfate precipitation
Concentration	Lot specific (U/mL)
Buffer	0.15 M phosphate buffered saline, pH 7.4
Preservative	0.01 % NaN <sub>3</sub>
Handling	Avoid freeze/thaw cycles
Long term storage	-18...-22°C
Expiration date	Specified on label