

PRODUCT DATA SHEET

Description	Monoclonal antibody to Alpha Fetoprotein
REF	10-01
Specificity	Reacts with human Alpha Fetoprotein
Clone	3G6, IgG1
Host Animal	Mouse
Immunogen	Human AFP
Cross-reactivity	Does not cross-react with serum albumin, chorionic gonadotropin, human placental lactogen and human immunoglobulin.
Source	Ascitic fluid
Format	Purified, liquid
Purification	Chromatography on Protein G Sepharose, purity is tested by electrophoresis
Concentration	lot-specific
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Long term storage	-18...-22°C
Handling	Avoid freeze/thaw cycles
Application	ELISA, use as detection antibody in ELISA with Mab 3H3