

PRODUCT DATA SHEET

Description	Monoclonal antibody to Alpha Fetoprotein
REF	10-02
Specificity	Reacts with human Alpha Fetoprotein
Clone	3H3, IgG1
Host Animal	Mouse
Immunogen	Human AFP
Cross-reactivity	Does not cross-react with serum albumin, chorionic gonadotropin, human placental lactogen and human immunoglobulin
Source	Ascitic fluid
Format	Purified, Liquid
Purification	Chromatography on Protein G Sepharose, purity is tested by electrophoresis
Concentration	Lot-specific
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium azide
Long term storage	-18...-22°C
Handling	Avoid freeze/thaw cycles
Application	ELISA, use as solid-phase antibody in ELISA together with Mab 3G6