

PRODUCT DATA SHEET

Description	Monoclonal antibody to Thyroxine
Specificity	Reacts with Thyroxine
Clone	3F3, IgG2b
Host Animal	Mouse
Immunogen	Thyroxine-BSA conjugate
Cross-reactivity, % L-Triiodothyronine	8.0
L-Diiodotyrosine	0.001
Source	Ascitic fluid
Format	Liquid
Purification	Not purified
Concentration	Lot-specific
Preservative	0.02% Sodium azide
Long term storage	-18...-22°C
Handling	Avoid freeze/thaw cycles
Application	ELISA
References	A Method for Cellular Stimulation of Thyroid Transplant for Correction of Thyroid System Function: Goldberg ED, Popov OS, Udut VV, Stavrova LA,Galyan AN, Zhdanov VV, Dygai AM <i>Bulletin of Experimental Biology and Medicine</i> , V. 139, № 4, 2005 , 529-531.

The human pituitary-thyroid system under chronic exposure to cold:Boiko ER, Potolitsyna NN, Kaneva AM

Doklady Biological Sciences, V. 416, № 1-6, 2007, 333-335.

On association between cortical 5-HT_{2A} receptors and behavior in rats with experimental thyroid disturbances:

Tikhonova MA, Kulikov AV, Lebedeva EI, Barykina NN, Amstislavskaya TG, Popova NK, *Pharmacology, biochemistry and behavior* V. 82, № 3, 2005, 506-514.

The pituitary-thyroid axis and oxygen consumption parameters under the conditions of chronic cold exposure in the North: Bojko EP, Evdokimov VG, Potolitsyna NN, Kaneva AM, Varlamova NG, Kochan TI, Vakhnina NA, Shadrina VD,

Solonin YG, Loginova TP, Eseva TV, Ketkina OA, Rogachevskaya OV, Lyudinina AY, Loginov AY

Human Physiology V. 34, № 2, 2008, 215-220.