

PRODUCT DATA SHEET

| | |
|--------------------|---|
| Description | Human IgE |
| REF | 11-04 |
| Source | Human myeloma cell line |
| Format | Liquid |
| Concentration | Lot specific (U/mL) |
| Buffer | Tris-acetate buffered saline, pH 7.4, contains glucose, glycerol, Pluronic-F68, BSA, Tartrazine |
| Preservative | 0.01 % NaN ₃ |
| Handling | Avoid freeze/thaw cycles |
| Long term storage | -18...-22 °C |
| Expiration date | Specified on label |